

RORA antibody (rAb), 100 µl

Catalog No: 91715

Clone: 230217G6

Application(s): WB

Reactivity: Human, Mouse, Rat

Volume: 100 µl

Purification: Protein A Chromatography

Host: Rabbit

Isotype: IgG

Molecular Weight: 59 kDa

Background: RORA is a steroid hormone receptor and has been implicated in the pathology of circadian rhythms, bone growth, angiogenesis, development of cones, cellular metabolism and a mediator in the immune and lipid metabolism pathways

Immunogen: This antibody was raised against the RORA fusion protein.

Buffer: Purified IgG in PBS, pH 7.3, with 50% glycerol and 0.02% sodium azide. Sodium azide is highly toxic.


Application Notes:

Applications Validated by Active Motif:

WB: 1:500-1:3000 dilution

Storage and Guarantee: Some products may be shipped at room temperature. This will not affect their stability or performance. Avoid repeated freeze/thaw cycles by aliquoting items into single-use fractions for storage at -20°C for up to 2 years. Keep all reagents on ice when not in storage. This product is guaranteed for 12 months from date of receipt.

This product is for research use only and is not for use in diagnostic procedures.



RORA antibody (rAb) tested by Western blot.

Detection of RORA by Western blot analysis using various lysates and probed with RORA antibody at 1:1500 dilution.